

Stress-Free Calf Tube Feeding



TRUSTI TUBER – STARTER PACK


UNIQUE DESIGN & EASY TO USE

Designed by a Vet, the Trusti Tuber significantly reduces calf stress and discomfort – soft, flexible tube, swallows easily and the calf struggling is minimised making the process very easy to learn.

- Flexible Tube can be held in place with one hand
- Large 4 Litre capacity bottle
- Transparent flexible tube makes it easy to clean & monitor flow
- Ergonomic bottle handle – perfect for large and small hands
- Reliable valve – fast flow even with thick colostrum
- Shown in trials to almost halve procedure feeding time.



Stopper & Calf size markers provide easy guide for tube placement & positioning



The mouthpiece sits comfortably against the muzzle.



Clip on bottle keeps mouthpiece clean & away from surfaces.

Flexible tube is advanced through mouthpiece so the tube is easily positioned within the oesophagus.



Comes with Peach Teat cap for easy transition between teat and tuber.

Contact us 1800 666 269 • daviesway.com.au

Product Code: 2500 – CALF TRUSTI TUBER STARTER PACK

DAVIESWAY
EST. 1911
IDEAS TO FARM. EVERYDAY.

HOW TO USE – TRUSTI TUBER



2 Pass the flexible tube through the mouthpiece so tube enters the oesophagus.



3 Stopper and calf size guides for easy positioning. When the tube has passed into the oesophagus it is safe to raise the bottle.



1 The mouthpiece sits comfortably against the muzzle

NUMBER ONE RULE:


Never let fluid flow unless the tube is safely in the oesophagus.

During Flow:

Hold the mouthpiece in place throughout feeding – hold with the calf so if the calf moves the mouthpiece moves with the calf. Do not let the tube pull out more than one inch during fluid flow. Never let the calf drop its head below the level of the stomach.


Removing the Tube:

Wait at least one second after all liquid has finished running through, then gently remove.



STOPPER POSITION & CALF SIZE MARKERS:

- Please consider the calf size markers a MINIMUM distance to pass the tube through.
 - The stopper should be positioned to the RIGHT HAND SIDE of the appropriate mark.
- Do not underestimate calf size as it is better to pass the tube farther than not far enough. If feeding a group of calves we recommend you set the stopper position to the largest calf in the mob.



Please watch the demonstration videos and read our calf handling tips – www.antahi.com



Daviesway Pty Ltd – 15 Helen Street, Heidelberg Heights, VIC 3081
Phone: +61 (0) 3 9459 4533 • 1800 666 269

For more information – daviesway.com.au

